

Micro-Plate Incubator

MPI2



Features

MPI2 is used for incubating the cell or enzymatic assay plates at appropriate temperature

- ➤ Micro-processing technology combined with PID
- > LCD display; temperature, time and speed display
- > Up and down, both side heating
- > Stable and reliable operation
- > Small, light and light noise
- > Alarm indication after stop operating
- ➤ Multiple safety protection, complying with CE, safe and convenient

Specifications	
Temp. Range	Ambient+5 ~ 80°C
Temp. Control Discrepancy	±0.3°C
Temp. Homogeneity	±0.3°C
Temp. Ascending Time	<15min (20° ~ 80°C)
Timer	0 ~ 99h59m
Capacity	2 micro-plates
Power	150W
Voltage	AC100~120V / AC200~240V; 50~60Hz
Weight	4Kg
Dimension (mm)	280(W) x 278(D) x 150(H)