

Dual-Gel Protein Electrophoresis

SPE1A



Features

SPE1A is a high quality polycarbonate molded SDS-PAGE unit for 2 gels; can run 1 gel or 2 gels at a time

- Exquisite molded bridge design
- Durable, no leakage
- High temperature resistance
- Accurate operation, save buffer
- High transparency, easy to observe
- Avoid bubble penetration between glass plates
- Detachable electrode holder and electrode head
- To ensure safety, power off while opening the cover
- Special Gel casting device for easy & fast gel casting
- Casting and running the gel at the same place
- Convenient to load samples and observe the bands
- Easy to thoroughly clean and repair

Specifications

Gel Size	83 x 75 mm; 2 Gels	
No. of Wells	10 & 15 wells	
Gel Thickness	1.0mm & 1.5mm	
Buffer Volume	400ml	
Dimension (mm)	197(L) x 96(W) x 100(H)	
Accessories Included	10 teeth comb (1.0mm) - 2 15 teeth comb (1.0mm) - 2 10 teeth comb (1.5mm) - 2 15 teeth comb (1.5mm) - 2	Glass plate set (1.0mm) - 4 Glass plate set (1.5mm) - 4 Blocking plate - 1 Casting Module - 1