

Dual-Gel Protein Electrophoresis

SPE1E



Features

SPE1E is a high quality SDS-PAGE unit for 2 gels of size 80 x 70mm each, for running up to 14 samples

- Durable, high temperature resistance
- Accurate operation, save buffer
- High transparency, easy to observe
- Convenient to load samples and observe the bands
- Withdrawable electrode, convenient for maintenance
- Designed for years of reproducible and rigorous use
- Easy and simple to use, save buffer solution
- Easy to thoroughly clean and repair
- Assortment of accessories supplied with the systems

Specifications

Gel Size	80 x 70 mm; 2 Gels	
No. of Wells	7 or 14	
Buffer Volume	150ml + 200ml	
Outer Dimension (mm)	160(L) x 200(W) x 150(H)	
Accessories Included	7 well Teflon comb (0.5mm) - 2 units 7 well Teflon comb (1.0mm) - 2 units Glass Plate (notched & Plain) - 2 sets Connecting Cord - 1 set (red & black)	Teflon Spacers (0.5mm) - 4 units Teflon Spacers (1.0mm) - 4 units Clamp & Screw - 2 sets Gasket - fixed