

Dual-Gel Protein Electrophoresis

SPE3



Features

SPE3 is a high quality polycarbonate molded SDS-PAGE unit for 2 gels; can run 1 gel or 2 gels at a time

- Exquisite molded bridge design
- Durable, no leakage
- High temperature resistance
- Accurate operation, save buffer
- High transparency, easy to observe
- Avoid bubble penetration between glass plates
- Detachable electrode holder and electrode head
- To ensure safety, power off while opening the cover
- Special Gel casting device for easy & fast gel casting
- Casting and running the gel at the same place
- Convenient to load samples and observe the bands
- Easy to thoroughly clean and repair

Specifications

Gel Size	130 x 100 mm; 2 Gels	
No. of Wells	12 & 16 wells	
Gel Thickness	1.0mm & 1.5mm	
Buffer Volume	1200ml	
Dimension (mm)	210(L) x 120(W) x 220(H)	
Accessories Included	12 teeth comb (1.0mm) - 2 16 teeth comb (1.0mm) - 2 12 teeth comb (1.5mm) - 2 16 teeth comb (1.5mm) - 2	Glass plate set - 4 Sealing Strip (1.0mm) - 4 Sealing Strip (1.5mm) - 4 Casting Module - 1