

## **Dual-Gel Protein Electrophoresis**

SPE4



## **Features**

SPE4 is a high quality polycarbonate molded SDS-PAGE unit for 2 gels; can run 1 gel or 2 gels at a time

- > Exquisite molded bridge design
- > Durable, no leakage
- > High temperature resistance
- > Accurate operation, save buffer
- ➤ High transparency, easy to observe
- > Avoid bubble penetration between glass plates
- > Detachable electrode holder and electrode head
- To ensure safety, power off while opening the cover
- > Special Gel casting device for easy & fast gel casting
- Casting and running the gel at the same place
- > Convenient to load samples and observe the bands
- Easy to thoroughly clean and repair

| Specifications       |  |  |
|----------------------|--|--|
| Gel Size             | 200 x 175 mm; 2 Gels   |  |
| No. of Wells         | 16 & 21 wells  |  |
| Gel Thickness        | 1.0mm  |  |
| Buffer Volume        | 3500ml   |  |
| Dimension (mm)       | 210(L) x 120(W) x 220(H)   |  |
| Accessories Included | 16 teeth comb (1.0mm) - 2<br>21 teeth comb (1.0mm) - 2<br>Glass plate set - 3<br>Sealing Strip (1.0mm) - 4 | Casting Module -1 Slanting Board -2 Silicon Tube - 2 PE Film - 1Pk |