

## **Tissue Homogenizer**

**TH15K** 



## **Features**

TH15K grinds the tissue samples in small centrifuge tubes and reduces sample loss during operation. It is a rechargeable device and supplied with 3.7 Volt, 600mAH Lithium battery. The battery can be fully charged in 2 hours through USB interface and last for over 5 hours.

- > Display is simple and clear
- Lightweight and easy to operate easy to use, energy saving and environmental protection
- Micro DC motor drives grinding pestle, which is strong power, short time and good effect
- Light and wear resistance, suitable for a variety of ways to eliminate the bacteria
- Quickly install the chuck, quick and easy disassembly
- > Supplied with 10 disposable grinding pestles, Optional stainless steel grinding pestle

| Specifications |               |
|----------------|---------------|
| Speed          | 0 ~ 15000rpm  |
| Noise          | ~ 50dB        |
| Dimension (mm) | 31(W) x 68(L) |
| Weight         | 200g          |