

Tissue Homogenizer

TH18K



Features

TH18K grind the tissue samples in small centrifuge tubes and reduces sample loss during operation

- Display is simple and clear
- Lightweight and easy to operate easy to use, energy saving and environmental protection
- Micro DC motor drives grinding pestle, which is strong power, short time and good effect
- Light and wear resistance, suitable for a variety of ways to eliminate the bacteria
- Quickly install the chuck, quick and easy disassembly
- Supplied with 10 disposable grinding pestles, Optional stainless steel grinding pestle

Specifications

Speed	8000 ~ 18000rpm
Noise	~ 60dB
Dimension (mm)	16(W) x 25(L)
Weight	400g