

Tissue Homogenizer

TH30K



Features

TH30K provides higher stability and good grinding performance. It is handy and easy to grind small amount of samples. The micro DC motor is used to drive grinding pestle for the wide applications. It's special for micro centrifuge tubes, and quick effect. Grinding pestle can be reused and probes are easily replaceable.

- Light and wear resistance, suitable for a variety of ways to eliminate the bacteria
- Lightweight and easy to operate easy to use, energy saving
- Quickly install the chuck, quick and easy disassembly
- The display is simple and clear
- Supplied with 10 disposable grinding pestles, Optional stainless steel grinding pestle

Specifications

Speed	8000 ~ 32000rpm
Noise	~ 72dB
Power	130W
Voltage	AC100~120V/AC200~240V; 50~60Hz
Weight	700g
Dimension (mm)	57(W) x 320(H)

Accessories

Model No.	Capacity	Stator Diameter	Rotor Diameter	Interval	Min. /Max. Depth	Length	Materials
TH30K/A	0.5-50ml	6mm	3.9mm	0.1mm	20/90mm	150mm	316L
TH30K/B	1-100ml	8mm	6.1mm	0.3mm	20/107mm	167mm	
TH30K/C	1-250ml	10mm	7.9mm	0.3mm	20/120mm	180mm	
TH30K/D	H-type stand; used to hold the homogenizer; height and angle are adjustable						